



Kolpen



AX

AX



AX.002.01/10
Dissecting Knives

AX.002.01 - AX.002.10

Dissecting Knives fig. 1-10
-All metal


Autopsy Knives, Cartilage Knives

SEZIERMESSER, KNORPELMESSER - CUCHILLOS PARA AUTOPSIA, CUCHILLOS PARA CARTÍLAGOS

Kolpen

GENERAL INSTRUMENTS
ALLGEMEINES INSTRUMENTARIUM
INSTRUMENTAL GENERAL

AX




AX.023.10 25,5 cm - 10"

1/2

AX.023.10

Autopsy Knives,
Cartilage Knives




1/2

AX.033.08


Virchow
Autopsy Knives,
Cartilage Knives

AX.023.10 | -Autopsy Knife, heavy pattern with wooden handle
AX.033.08 | -Cartilage Knife

AX



AX.036.11/17
Walz
Post Mortem Knives,
Brain Knives



AX.038.16
AX.038.20
AX.038.24



AX.038.16/24
Virchow
Post Mortem Knives,
Brain Knives

AX.036.11 - AX.036.17

AX.038.16 - AX.038.24

-Post Mortem Knives (Organ Knives), with
wooden handles
-Brain Knives, with metal handles

Razor Blade Holders (and Breakers), Razor Blade, Micro Blade

RASIERKLINGENHALTER UND- BRECHER, RASIERKLINGEN, MIKROKLINGEN
 PINZAS PARA SUJETAR Y ROMPER HOJAS DE AFEITAR, HOJA DE AFEITAR,
 MICROHOJA

Kolpen

GENERAL INSTRUMENTS
 ALLGEMEINES INSTRUMENTARIUM
 INSTRUMENTAL GENERAL


AX



AX.090.09
 10,0 cm - 4"
 Ø 5,5 mm

1/1

AX.090.09
 Troutman-Chris
 Razor Blade Holders



AX.091.09
 10,0 cm - 4"
 Ø 5,5 mm

1/1


AX.091.09
 Troutman-Chris
 Razor Blade Holders



AX.092.11
 11,0 cm - 4 1/2"

1/1

AX.092.11
 Razor Blade Holders



AX.094.13
 13,0 cm - 5 1/4"
 Ø 8 mm

1/1

AX.094.13
 Castroviejo
 Razor Blade Holders



AX.096.00
 Package of 10 Blades
 Razor Blades



AX.097.00
 Individually packed,
 sterile, in stripes of 10
 Micro Blades


1/1

- AX.096.00 | -Razor Blade, breakable, stainless steel
- AX.097.00 | -Micro Blade, 4/100 mm thick, specially suited for cataract surgery

Also available of Titanium

AX

AX.100.04
Solid



1/1

AX.100.04
Standard No. 4
Scalpel Handles

AX.100.34
Double
Handle No. 4
+ No. 3


Solid



1/1


AX.100.34
Scalpel Handles

AX.104.41
Solid
Straight



1/1

AX.104.41
No. 4 L
Scalpel Handles



AX.104.45
Solid
Angled

1/1

AX.104.45
No. 4 L
Scalpel Handles

Scalpel Handles

SKALPELL-GRIFFE - MANGOS DE BISTURÍES

Kolpen


GENERAL INSTRUMENTS
ALLGEMEINES INSTRUMENTARIUM
INSTRUMENTAL GENERAL

AX

For Blades

AX.150.18 - AX.151.36


AX.155.18 - AX.156.36



AX.108.08
Hollow

AX.108.08

Collin
No. 8
Scalpel Handles



AX.109.09
Hollow

AX.109.09

Collin
No. 9
Scalpel Handles



AX.110.10
Hollow

AX.110.10

Martin-Super
No. 10
Scalpel Handles



AX.112.04
Plastics


AX.112.04

Scalpel Handles

AX

For Blades
AX.150.10 - AX.150.17
AX.151.12
AX.151.15
AX.155.09
AX.155.10 - AX.155.17
AX.156.12
AX.156.15

AX.133.31
Solid
Straight




AX.133.35
Solid
Angled



AX.130.03
Solid

1/1

AX.130.03
Standard No. 3
Scalpel Handles



AX.130.04
Solid

1/1

AX.130.04
No. 3
Scalpel Handles



AX.130.05
Solid


1/1

AX.130.05
No. 3
Scalpel Handles



1/1

AX.133.31
No. 3 L
Scalpel Handles



1/1

AX.133.35
No. 3 L
Scalpel Handles

Scalpel Handles

SKALPELL-GRIFFE - MANGOS DE BISTURÍES

Kolpen

GENERAL INSTRUMENTS
ALLGEMEINES INSTRUMENTARIUM
INSTRUMENTAL GENERAL

AX



AX.137.07
Solid

1/1

AX.137.07
No. 7
Scalpel Handles



AX.137.71
Solid

1/1

AX.137.71
No. 7 K
Scalpel Handles



AX.140.03
Hollow,
round

1/1

AX.140.03
Scalpel Handles





AX.142.03
Plastics

1/1

AX.142.03
Scalpel Handles

AX



Non-Sterile

Sterile

Saw Blade



AX.155.09

Single Sterile Package
Scalpel Blades




Fig. 10

AX.150.10

Fig. 10

AX.155.10



Fig. 11

AX.150.11

Fig. 11

AX.155.11



Fig. 12

AX.150.12

Fig. 12

AX.155.12



Fig. 13

AX.150.13

Fig. 13

AX.155.13



Fig. 14

AX.150.14

Fig. 14

AX.155.14



Fig. 15

AX.150.15

Fig. 15

AX.155.15

Scalpel Blades


SKALPELL-KLINGEN - HOJAS DE BISTURÍES

Kolpen

GENERAL INSTRUMENTS
ALLGEMEINES INSTRUMENTARIUM
INSTRUMENTAL GENERAL

AX

- For Handles
- AX.100.34
 - AX.130.03
 - AX.130.04
 - AX.133.31
 - AX.133.35
 - AX.137.07
 - AX.137.71
 - AX.140.03
 - AX.142.03



* Double cutting edge




	Fig. 16	Fig. 17	Fig. 12d	Fig. 15c
Non-Sterile	AX.150.16	AX.150.17	AX.151.12	AX.151.15
Sterile	AX.155.16	AX.155.17	AX.156.12	AX.156.15

AX

For Handles AX.100.04 AX.100.34 AX.104.41 AX.104.45
AX.108.08 AX.109.09 AX.110.10 AX.112.04







Fig. 18 Fig. 19 Fig. 20 Fig. 21 Fig. 22 Fig. 23 Fig. 24

Non-Sterile

AX.150.18 AX.150.19 AX.150.20 AX.150.21 AX.150.22 AX.150.23 AX.150.24

Sterile

Fig. 18 Fig. 19 Fig. 20 Fig. 21 Fig. 22 Fig. 23 Fig. 24

AX.155.18 AX.155.19 AX.155.20 AX.155.21 AX.155.22 AX.155.23 AX.155.24

Scalpel-Blade Remover



AX.199.00

Scalpel Blades

* Double cutting edge





	Fig. 25	Fig. 34	Fig. 36	Fig. 24d	Fig. 36d
Non-Sterile	AX.150.25	AX.150.34	AX.150.36	AX.151.24	AX.151.36
Sterile	AX.155.25	AX.155.34	AX.155.36	AX.156.24	AX.156.36

AX

Sterile



1/1

AX.285.18

Fig. 18
Disposable
Scalpels
with

Plastic Handle

AX.285.19

Fig. 19
Disposable
Scalpels
with

Plastic Handle

AX.285.20

Fig. 20
Disposable
Scalpels
with

Plastic Handle

AX.285.21

Fig. 21
Disposable
Scalpels
with

Plastic Handle

AX.285.22

Fig. 22
Disposable
Scalpels
with

Plastic Handle

AX.285.23

Fig. 23
Disposable
Scalpels
with

Plastic Handle

AX.285.24

Fig. 24
Disposable
Scalpels
with

Plastic Handle

AX.285.25

Fig. 25
Disposable
Scalpels
with


Plastic Handle

Disposable Scalpels with Plastic Handle

EINMAL-SKALPELLE MIT KUNSTSTOFF-GRIFF - BISTURÍES PARA USO ÚNICO
 CON MANGOS DE MATERIA PLÁSTICA

Sterile

* Double cutting edge



AX.295.10
 Fig. 10
 Disposable
 Scalpels
 with
 Plastic Handle

AX.295.11
 Fig. 11
 Disposable
 Scalpels
 with
 Plastic Handle

AX.295.12
 Fig. 12
 Disposable
 Scalpels
 with
 Plastic Handle

AX.295.13
 Fig. 13
 Disposable
 Scalpels
 with
 Plastic Handle

AX.295.15
 Fig. 15
 Disposable
 Scalpels
 with
 Plastic Handle

AX.296.12
 Fig. 12d
 Disposable
 Scalpels
 with
 Plastic Handle

AX.296.15
 Fig. 15c
 Disposable
 Scalpels
 with
 Plastic Handle

AX

AX.410.01
Solid



1/1

AX.410.01
Operating Knives,
Amputating Knives

AX.500.13 - AX.510.17 | -Amputating Knives
AX.520.12 | -Flap Knife
AX.530.09 | -Phalangeal Knife
AX.540.11 | -Interosseous Knife



1/2

Lengths of Blades				
cm	in.	Catlin	Liston	Langenbeck
9,5	3 3/4			AX.530.09
11,0	4 1/2			AX.540.11
12,00	4 3/4			AX.520.12
13,00	5	AX.500.13		
16,0	6 1/4	AX.500.16	AX.506.16	
17,0	6 3/4			AX.510.17
19,0	7 1/2	AX.500.19	AX.506.19	
22,0	8 3/4	AX.500.22	AX.506.22	
28,0	11 1/4	AX.500.28		

Resection Knives, Razors

RESEKTIONSMESSER, RASIERMESSER -
CUCHILLOS PARA RESECCIÓN, NAVAJAS

Kolpen

GENERAL INSTRUMENTS
ALLGEMEINES INSTRUMENTARIUM
INSTRUMENTAL GENERAL

AX




AX.570.05
Cutting Length:
5,5 cm - 2 1/4"

1/2

AX.570.05

Langenbeck
Resection Knives,
Razors



AX.572.05
Cutting Length:
5,5 cm - 2 1/4"

1/2

AX.572.05

Langenbeck
Resection Knives,
Razors



AX.800.00
Cutting Length:
7,5 cm - 3"

1/1

1/2

AX.800.00

Resection Knives,
Razors



AX.811.00
Cutting Length:
8,0 cm - 3 1/4"

1/1

1/2

AX.811.00

Resection Knives,
Razors



AX.870.00
210 x 30 mm

AX.870.00


Resection Knives,
Razors

AX.890.00
100 mm

AX.890.00

Resection Knives,
Razors

AX.875.00
150 x 50 mm



AX.875.00

Record-Arkansas
Resection Knives,
Razors

AX.800.00
AX.870.00


-Both sides hollow ground, metal handle
-Sharpening Strop, 4 sides, one side finest emery stone, one side prepared leather, 2 sides red and black prepared

AX.875.00
AX.890.00

-Arkansas Oil Stone, in wooden box
-Silicon-Carbide Grindstone

AX

AX.880.00
30,0 cm - 12"



1/2


AX.880.00
Schink
Dermatomes

1/2




AX.880.01
Dermatomes

1/2



AX.880.02
Dermatomes



1/2

AX.882.19
Silver
Dermatomes

AX.880.00


AX.880.01

AX.880.02

-Dermatome, adjustable from 0,1 to 2,0 mm, compl., with 1 interchangeable blade and 2 skin straightening blades

-Spare Blade only

-Skin straightening plate only



AX.900.04 16,0 cm - 6 1/4"
cutting diameter 4 mm

AX.900.05 16,0 cm - 6 1/4"
cutting diameter 5 mm

AX.900.04/05
Aorta Punches